Luminaire for Diagnosis

LUVIS | EI100

* Be sure to read the manual before using this product.

1 Standard	03	2 Cautions and Warnings	04
3 Symbol	08	4 Introduction	10
5 Product specification	14	4.1 Intended use 4.2 General description 4.3 Environmental requirement 4.4 Safety information 4.5 Warranty regulation 4.6 Electromagnetic Compatibility	
5.1 Technical specification (In accordance with EN/IEC 60601-2-41) 5.2 Electrical specification (In accordance with EN/IEC 60601-1) 5.3 Mechanical specification		6 Other characteristic	14
7 Use	15	7.1 HEAD CONTROLLER 7.2 Positioning	
		8 Product installation	16
		8.1 HEAD and MAIN ARM 8.2 E100 Mobile(E1-M) 8.3 E100 Clamp(E1-C) 8.4 E100 Wall(E1-W) 8.5 Checking cover and cap 8.6 Operating range	
9 Cleaning	19	10 Maintenance	19
11 Troubleshooting	20	12 Disposal	20
13 Model designation	21	14 List of components	21
15 Manufacturer and factory	22	16 Europe representative	23

1. Standard

* Certification of DENTIS

- EN ISO 13485:2016
- Relevant EC Regulation: REGULATION (EU) 2017/745

* Applied Standards

- EN ISO 15223-1:2016, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
- EN 1041:2008/A1:2013, Information supplied by the manufacturer with medical devices
- EN ISO 13485:2016, Medical devices – Quality management systems – Requirements for regulatory purpose
- EN ISO 14971:2019, Medical devices – Application of risk management to medical devices
- EN 60601-1:2006/A1:2013, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
- EN 60601-1-2:2015, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances – Requirements and tests
- EN 60601-1-6:2010, Medical electrical equipment-Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability
- EN ISO 7010:2019, Graphical symbols – Safety colors and safety signs-Registered safety signs
- EN 62471:2008, Photo biological safety of lamps and lamps systems
- EN 62366-1:2015, Medical devices – Application of usability engineering to medical devices
- IEC 60601-2-41:2009/A1:2015, Medical electrical equipment – Part 2-41: Particular requirements for the basic safety and essential performance of surgical luminaires and luminaires for diagnosis
- IEC 62133:2012, Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications
- IEC 62304:2006/A1:2015, Medical device software, Software life-cycle processes

2. Cautions and Warnings

⚠ CAUTION(ATTENTION) This ME Equipment is intended only for use in the professional health-care facility environment, e.g. public and private hospitals, specialized medical offices, etc. This ME Equipment is intended for direct use on operation theatres and nearby HF surgical equipments, where the intensity of EM DISTURBANCES complies with the applicable standards. For further instructions please follow the chapter "Electromagnetic compatibility" of the User's Manual.

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⚠ CAUTION(ATTENTION) If this ME Equipment is lost or degraded the performance due to EM DISTURBANCES, result in improper operation and degradation of the performance.

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⚠ CAUTION(ATTENTION) All cables and maximum lengths of coaxial cables that are replaceable by the DENTIS and that are likely to affect compliance of the this ME Equipment with the requirements of EMC(Electro-Magnetic compatibility). Do not modify this ME Equipment.

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⚠ CAUTION(ATTENTION) The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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⚠ WARNING(AVERTISSEMENT) Use of this ME Equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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⚠ WARNING(AVERTISSEMENT) Use of accessories, transducers and cables other than those specified or provided by DENTIS of this ME Equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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⚠ WARNING(AVERTISSEMENT) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm(12 inches) to any part of the E100, including cables specified by DENTIS. Otherwise, degradation of the performance of this equipment could result.

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⚠ WARNING(AVERTISSEMENT) The instructions given in this document must be followed when handling the product. Failure to do so may endanger the safety of the installers or users.

As well as specific information on operating the entire product and conducting preventive maintenance, are provided in the USER'S MANUAL.

For further information, please contact our sales network or our local network.

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⚠ WARNING(AVERTISSEMENT) The electrical connections must be performed by a qualified technician only.

The electrical installation must be planned, performed and inspected by electrical engineers.

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⚠ WARNING(AVERTISSEMENT) The LIGHTHEAD is designed to operate using a AC 100-240V, 50-60 Hz. Higher or lower voltages may affect the light intensity and operating life of the LEDs. Damaged wire insulation may result in the risk of electric shock.

To protect against the risk of electrocution, connect the power cables carefully.

Luminaire for diagnosis is class II for protection against electric shock.

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⚠ WARNING(AVERTISSEMENT) Damaged wire insulation may result in the risk of electric shock.

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⚠ WARNING(AVERTISSEMENT) This product may only be repaired and special assembly work may only be carried out by DENTIS or a company that has been authorized by DENTIS.

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⚠ WARNING(AVERTISSEMENT) Do not look directly into light source(LED).

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⚠ WARNING(AVERTISSEMENT) All the information in this manual has been checked out carefully and discerned as accurate one at the time of publication.

However, DENTIS takes no responsibilities of the results caused by default, omission, or misuse of it.

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⚠ WARNING(AVERTISSEMENT) DENTIS has rights to modify the product itself or specifications of the product without any prior notice, as well as rights not to renew that modification on this manual.

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⚠ WARNING(AVERTISSEMENT) Do not press more than two buttons simultaneously.

In case of abnormal operation (overpower) of this product, stop the medical treatment and contact the place of purchase.

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
⚠ WARNING(AVERTISSEMENT) The circuit of medical device must be installed in the state with the means to electrically separate with in all poles from the POWER SUPPLY.

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⚠ WARNING(AVERTISSEMENT) When using the Adapter type, be sure to make easy connect and disconnect the power cord.

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No	Location
①	Marking on the outside of ME EQUIPMENT
②	Marking on the inside of ME EQUIPMENT
③	Marking on the controls of ME EQUIPMENT
④	Marking on the packing label of ME EQUIPMENT
⑤	Marking on the manual of ME EQUIPMENT
⑥	Marking on the label of ME EQUIPMENT

4. Introduction

4.1 Intended use

- Luvis Luminaire for Diagnosis(E100) is used for treatment or diagnosis in operating theatres or in medical rooms (groups 0, 1 and 2 according to IEC 60364-7-710)
- Classification under the provision of 93/42/EEC(MDD) : Class I
 - Luvis Luminaire for Diagnosis(E100) is classified as a Class I device.
- Form of protection against electric shock : Class II
 - Luvis Luminaire for Diagnosis(E100) is classified as Class II
- Degree of protection against flammability
 - Luvis Luminaire for Diagnosis(E100) is classified as a device not suitable to be used in a potentially flammable environment.
 - Do not use near flammable materials.
- Method(s) of cleaning recommended by the manufacturer.
 - The LIGHTHEAD should be cleaned with cloths and a specified disinfectant.
- Mode of operation
 - Classification of Luvis Luminaire for Diagnosis(E100) : continuous operation.

4.2 General description

- The user must ensure that the device works properly and is in a satisfactory condition before each use.
- This DENTIS device is intended only for use in the field of medical. It is impermissible to use the device for a purpose for which it was not intended.
- "Proper Use" includes following all the instructions for use and ensuring that all inspections and service tasks are performed.
- Apply and meet the overarching guidelines and/or national laws, national regulations and the rules of technology for medical devices applicable for startup and use of the device for the intended purpose.
- The user must observe the following:
 - Only use properly operating equipment.
 - Protect himself or herself and third parties from danger.
 - Avoid contamination from the device.
- During use, the following national regulations must be observed:
 - The applicable health and safety regulations.
 - The applicable accident prevention regulations.
- To ensure that device maintains their value and are always ready for use, they must be serviced once a year as recommended.
- The safety checks must be performed every year.
- Repair and service of the device is authorized only to those who meet the requirements below:
 - Technicians of authorized dealers specially trained by DENTIS.
 - The trained technicians of DENTIS branches.

4.3 Environmental requirement

- Conditions of the usage environment
 - Temperature : 0 - 40°C
 - Relative Humidity : 30 - 90%RH
 - Atmospheric pressure : 0-2,000 m (106-80 kpa)
- Conditions of storage and transportation environment
 - Temperature : 0 - 40 °C
 - Relative Humidity : 0 - 80 %
 - Atmospheric pressure: 0 - 2,000 m (106 – 80 Kpa)

4.4 Safety information

- Luvis Luminaire for Diagnosis(E100) as a medical device complies with the safety regulation EN/IEC 60601-1.
- Furthermore all configurations shall comply with the system standard EN 60601-1-2:2015.
- If in doubt, consult the technical service department or your local representative.
- For EU Countries
Europe Representative
KTR Europe GmbH
Mergenthalerallee 77, 65760 Eschborn, Germany
- Safety is everyone's obligation and responsibility.
- The safe use of this product is related to everyone such as installer, user, operator and equipment's manager.
- It must read and learn this user's manual is compulsory before installation, using, cleaning, fixing or operation of this product or its accessories. Pay particular attention and be familiar with warning symbols about safety.
- If do not follow safety direction of this manual, you can get injured or accident when you operate this product. After read carefully and understand this manual, use this product.
- This manual is in keep a place where you can find easily.

4.5 Warranty regulation

- DENTIS warrants all products against defects in materials or workmanship for two year from time of delivery.
- DENTIS's sole obligation under the warranty is to provide parts for the repair, or at its option, to provide the replacement product (excluding labor).
- The buyer shall have no other remedy. All special, incidental, and coincidental damages are excluded.
- Written notice of breach must be given to DENTIS within the warranty period.
- The warranty does not cover damage resulting from improper installation or maintenance, accident or misuse.
- The warranty does not cover damage resulting from the use of cleaning, disinfecting or sterilization chemicals and processes.
- The Failure to follow instructions provided in the DENTIS Instructions for Use (operation and maintenance instructions) may void the warranty.
- LED PCB Ass'y is covered under 60,000 hrs warranty.

4.6 Electromagnetic Compatibility

4.6.1 Emission

This ME equipment is intended for use in Professional healthcare facility environment.		
Emission test	Compliance	Guidance
Conducted Disturbance CISPR 11(EN 55011)	Complies (Group 1, Class A)	This ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated Disturbance CISPR 11(EN 55011)		
Harmonic current IEC 61000-3-2	Complies	This ME equipment is suitable for use in all establishments other than domestic premises and those directly connected to the public low-voltagepower supply network that supplies buildings used for domestic purposes.
Voltage fluctuations and flicker IEC 61000-3-3	Complies	

4.6.2 Immunity

This ME equipment is intended for use in Professional healthcare facility environment.																										
Immunity test	EN 60601-1-2:2015	Compliance																								
Electrostatic Discharge(ESD) IEC 61000-4-2	Direct: ± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air Indirect: ± 8 kV HCP/VCP	Complies																								
Radio Frequency Electromagnetic Fields IEC 61000-4-3	3 V/m @ 80 MHz ~ 2.7 GHz 80 % AM at 1 kHz	Complies																								
Proximity fields from RF wireless communications equipment IEC 61000-4-3	<table border="1"> <thead> <tr> <th>Frequency (MHz)</th> <th>Modulation</th> <th>Immunity Level (V/m)</th> </tr> </thead> <tbody> <tr> <td>385</td> <td>**Pulse Modulation: 18 Hz</td> <td>27</td> </tr> <tr> <td>450</td> <td>*FM + 5Hz deviation: 1 kHz sine</td> <td>28</td> </tr> <tr> <td>710 745 780</td> <td>**Pulse Modulation: 217 Hz</td> <td>9</td> </tr> <tr> <td>810 870 930</td> <td>**Pulse Modulation: 18 Hz</td> <td>28</td> </tr> <tr> <td>1 720 1 845 1 970</td> <td>**Pulse Modulation: 217 Hz</td> <td>28</td> </tr> <tr> <td>2 450</td> <td>**Pulse Modulation: 217 Hz</td> <td>28</td> </tr> <tr> <td>5 240 5 500 5 785</td> <td>**Pulse Modulation: 217 Hz</td> <td>9</td> </tr> </tbody> </table>	Frequency (MHz)	Modulation	Immunity Level (V/m)	385	**Pulse Modulation: 18 Hz	27	450	*FM + 5Hz deviation: 1 kHz sine	28	710 745 780	**Pulse Modulation: 217 Hz	9	810 870 930	**Pulse Modulation: 18 Hz	28	1 720 1 845 1 970	**Pulse Modulation: 217 Hz	28	2 450	**Pulse Modulation: 217 Hz	28	5 240 5 500 5 785	**Pulse Modulation: 217 Hz	9	Complies
	Frequency (MHz)	Modulation	Immunity Level (V/m)																							
	385	**Pulse Modulation: 18 Hz	27																							
	450	*FM + 5Hz deviation: 1 kHz sine	28																							
	710 745 780	**Pulse Modulation: 217 Hz	9																							
	810 870 930	**Pulse Modulation: 18 Hz	28																							
	1 720 1 845 1 970	**Pulse Modulation: 217 Hz	28																							
	2 450	**Pulse Modulation: 217 Hz	28																							
	5 240 5 500 5 785	**Pulse Modulation: 217 Hz	9																							
	** The carrier shall be modulated using a 50 % duty cycle square wave signal. * As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.																									
Fast Transients IEC 61000-4-4	<table border="1"> <thead> <tr> <th>Voltage</th> <th>AC/DC power ports</th> <th>Signal ports</th> </tr> </thead> <tbody> <tr> <td>Test voltage</td> <td>± 2 kV</td> <td>± 1 kV</td> </tr> </tbody> </table>	Voltage	AC/DC power ports	Signal ports	Test voltage	± 2 kV	± 1 kV	Complies																		
Voltage	AC/DC power ports	Signal ports																								
Test voltage	± 2 kV	± 1 kV																								
- 100 kHz repetition frequency																										
Surges IEC 61000-4-5	<table border="1"> <thead> <tr> <th>Voltage</th> <th>Power lines</th> </tr> </thead> <tbody> <tr> <td>Test voltage</td> <td>Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV</td> </tr> </tbody> </table>	Voltage	Power lines	Test voltage	Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV	Complies																				
Voltage	Power lines																									
Test voltage	Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV																									
RF Continuous Conducted IEC 61000-4-6	3 V @ 0.15 MHz ~ 80 MHz 6 V @ in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Complies																								
Power Frequency Magnetic Fields IEC 61000-4-8	30 A/m @ 50 Hz or 60 Hz	Complies																								
Voltage Dips, Interruptions, and Variations IEC 61000-4-11	<ul style="list-style-type: none"> • Voltage Dips 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° • Voltage Interruptions 0 % UT; 250/300 cycle • Voltage Variations 	Complies																								
	<table border="1"> <thead> <tr> <th>Frequency (Hz)</th> <th>Ranges</th> </tr> </thead> <tbody> <tr> <td>50</td> <td>49, 50, 51</td> </tr> <tr> <td>60</td> <td>59, 60, 61</td> </tr> </tbody> </table>		Frequency (Hz)	Ranges	50	49, 50, 51	60	59, 60, 61																		
	Frequency (Hz)		Ranges																							
50	49, 50, 51																									
60	59, 60, 61																									

5. Product Specification

5.1 Technical specification (In accordance with EN/IEC 60601-2-41)

· STANDARD SPECIFICATION

Specifications	Unit	E100	Remark
Central illuminance Ec (@50cm)	lx	30,000 / 45,000 / 60,000	
Electronic adjustment range	%	50~100	
Light field diameter d10	cm	17	
colour temperature (3 levels)	K	3,500 / 4,000 / 4,500	
colour rendering index, Ra	N/A	95	

* Optical values are measured with a tolerance of ±10%

5.2 Electrical specification (In accordance with EN/IEC 60601-1)

· ADAPTER TECHNICAL DATA

Content	Luvis E100
Input	AC 100~240 V, 50-60 Hz, 0.6 A
Output	DC 15 V, 1.67 A
Power	Max 25 W

Wall Plug-in SMPS
(Interchangeable Blades)

· E100 RATING

Content	Rating	Remark
E100	AC 100-240 V, 28-46 VA, 50-60 Hz	

⚠ WARNING : Use the special adapter, which Must provided by the DENTIS.

5.3 Mechanical specification

· Component

Content	Weight(kg)
LIGHTHEAD + FIRST HEAD ARM + SECOND HEAD ARM	1.9
Mobile Stand (VERTICAL ARM + BASE)	4.9
CLAMP	0.5
WALL BRACKET	0.5
CEILING BRACKET	1.0
CEILING ARM	4.6


6. Other characteristic

Specifications	Luvis E100	Remark
Protection against electric shock	Class II Protection	
Protection against harmful ingress of water or particulate matter	Ordinary	Head part : IP43

7. Use

7.1 HEAD CONTROLLER


· The LIGHT HEAD can be controlled using the HEAD CONTROLLER

No.	Name	Description
①	Standby Button	Power On/Off
②	Intensity Control Button	3 Steps adjustable : 30,000/45,000/60,000 lx
③	Color Temp. Control Button	3 Steps adjustable : 3,500/4,000/4,500 K

7.2 Positioning


· Positioning the E100 LIGHTHEAD

8. Performance characteristic

8.1 HEAD and MAIN ARM


· HEAD and ARM PART Tension setting

* Loosen the KNOB to adjust the position and tighten the KNOB to fix the position when moving it up and down


8.2 E100 Mobile(E1-M)

· VERTICAL ARM Height setting

* Loosen the KNOB to adjust the position and tighten the KNOB to fix the position when moving it up and down

8.3 E100 Clamp(E1-C)


* Loosen the KNOB to adjust the position and tighten to fix the position.

8.4 E100 Wall(E1-W)

* Loosen the KNOB to adjust the position and tighten the KNOB to fix the position.

8.5 Checking cover and cap



* Check that all covers and bolt caps are in place.

8.6 Operating range of E100


· MOBILE TYPE

(Unit : mm)


· CLAMP TYPE

(Unit : mm)


· WALL TYPE

(Unit : mm)

· CEILING TYPE

(Unit : mm)

9. Cleaning

⚠ WARNING

Before and after use, clean the LIGHTHEAD with alcohol or ethanol.
It is not allowed to use agents which contain the following substances :
High-concentration organic and inorganic acids, Chlorinated hydrocarbons.

10. Maintenance


⚠ WARNING

1. A electrical and mechanical check-up should be done every year.
Disconnect the pendant system from the mains before any maintenance work to prevent electric shock.
2. Check points on maintenance work
 - Defects of paint
 - Plastic fissures
 - Loosened parts.
 - Free rotatability / limit stops and deformation of the suspension
 - The connection between LIGHTHEAD and components
 - Defect of the KNOB
3. In case of failure or damage, please contact your supplier.

⚠ WARNING

1. Do not place things on or hang on MAIN ARM.
2. If this device is operated beyond the indicated angles, it can be damaged.


11. Troubleshooting

No.	Problem	Cause	Corrective action
1	the LED light on the LIGHTHEAD does not turn on.	Power cut	Check if supply mains are operating.
		Others	Press the stand-by button of the control panel and check that the LED turns on. (If green light does not light on, do not attempt to repair it yourself, please contact your supplier)
2	If LIGHTHEAD is not fixed after handling, moves up and down.	KNOB incorrectly adjusted Loosen the KNOB to adjust the position and tighten the KNOB to fix the position when moving it up and down	 <p>As shown in the picture, Immerse the KNOB tightly to fix the position</p>
3	If there is an error in the light field diameter.	Inappropriate distance	Check if recommended distance(50cm) between the LIGHTHEAD to the patient's body. If the condition does not improve, do not attempt to repair. Please contact your supplier.
4	If the SAFETY COVER of the LIGHTHEAD is polluted.	Pollution	Cleanse it with designated chemicals (Alcohol, Ethanol). If the essential performances of the device (Intensity of Illumination, colour temperature) are seriously impaired, please contact your supplier.

13. Model designation

Model	-	Mount Type
E100	-	E=Ceiling Type W=Wall Type C=Clamp Type M=Mobile Type

14. List of components


LIGHTHEAD + HEAD ARM	MOBILE VERTICAL ARM	CLAMP	WALL BRACKET	CEILING BRACKET & ARM HOUSING CAP
				
CEILING VERTICAL ARM	CEILING COVER and BRACKET	ADAPTER		
				

12. Disposal

-  1. For environment and safety of human, wastes must be recycled or separated.
- The materials should be carefully separated.
 - The electrical boards should be submitted to an appropriate recycling proceeding
 - The cardboard box may be recycled with other paper products.
2. Please contact the closest DENTIS branch or your supplier, if you have any questions about recycling of the device.


15. Manufacturer and Factory

Manufacturer, Factory : DENTIS MEDICAL DIVISION
6, Yuram-ro, Dong-gu, Daegu, Republic of Korea

16. Europe Representative

Location of KTR Europe GmbH. Mergenthalerallee 77, 65760 Eschborn, Germany

DENTIS

DENTIS Co., Ltd. (HEADQUARTERS) 87, Yuram-ro Dong-gu, Daegu, Republic of KOREA Tel +82-53-583-2804 Fax +82-53-583-2806
DENTIS MEDICAL DIVISION(FACTORY) 6, Yuram-ro, Dong-gu, Daegu, Korea Tel. +82-70-7730-4146 Email. luvis@luvis.co.kr

EL10GP-PM030