Minor Surgical Luminaire

LUVIS | S200

* Be sure to read the manual before using this product.

1 Standards	03	2 Cautions and Warnings	04
3 Symbol	08	4 Introduction	09
5 LIGHTHEAD specification	13	4.1 Intended use	
5.1 Product colors : BLUE, IVORY, WHITE		4.2 General description	
5.2 Technical specification (In accordance with EN/IEC 60601-2-41)		4.3 Environmental requirement	
5.3 Electrical specification (In accordance with EN/IEC 60601-1)		4.4 Safety information	
5.4 Mechanical specification		4.5 Warranty regulation	
5.5 Other characteristic		4.6 Electromagnetic Compatibility	
		4.6.1 Emission	
		4.6.2 Immunity	
7 Labelling requirement	20	6 USE	15
7.1 MARKING PLATE		6.1 Product type	
7.2 CAUTION AND INFORMATION MARK		6.2 Product description	
		6.3 Performance characteristic	
		6.4 Detachable part	
8 Sterilization and Cleaning	22	9 Operation	25
8.1 Sterilization			
8.2 Cleaning			
10 Adjustment and Maintenance	28	11 Disposal	29
10.1 Adjustment			
10.2 Maintenance			
13 Manufacturer and Factory	31	12 Model designation	30
		14 Europe Representative	31

1. Standards

* Certification of DENTIS

- EN ISO 13485:2016
- Relevant EC Regulation: REGULATION (EU) 2017/745

* Applied Standards:

- EN ISO 15223-1:2016, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.
- EN 1041:2008/A1:2013, Information supplied by the manufacturer with medical devices.
- EN ISO 13485:2016, Medical devices – Quality management systems – Requirements for regulatory purpose.
- EN ISO 14971:2019, Medical devices – Application of risk management to medical devices.
- EN 60601-1:2006+A2:2021, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.
- EN 60601-1-2:2015, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances – Requirements and tests.
- EN 60601-1-6:2010, Medical electrical equipment-Part 1-6: General requirements for basic safety and essential performance Collateral Standard: Usability.
- EN ISO 7010:2019, Graphical symbols – Safety colors and safety signs-Registered safety signs.
- EN 62471:2008, Photo biological safety of lamps and lamps systems.
- EN 62366-1:2015, Medical devices – Application of usability engineering to medical devices.
- IEC 60601-2-41:2009/A1:2015, Medical electrical equipment – Part 2-41: Particular requirements for the basic safety and essential performance of surgical luminaires and luminaires for diagnosis.
- IEC 62304:2006/A1:2015, Medical device software, Software life-cycle processes.

2. Cautions and Warnings

⚠ CAUTION(ATTENTION) This ME Equipment is intended only for use in the professional health-care facility environment, e.g. public and private hospitals, specialized medical offices, etc. This ME Equipment is intended for direct use on operation theatres and nearby HF surgical equipments, where the intensity of EM DISTURBANCES complies with the applicable standards. For further instructions please follow the chapter "Electromagnetic compatibility" of the User's Manual.

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⚠ CAUTION(ATTENTION) If this ME Equipment is lost or degraded the performance due to EM DISTURBANCES, result in improper operation and degradation of the performance.

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⚠ CAUTION(ATTENTION) All cables and maximum lengths of coaxial cables that are replaceable by the DENTIS and that are likely to affect compliance of the this ME Equipment with the requirements of EMC(Electro-Magnetic compatibility). Do not modify this ME Equipment.

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⚠ CAUTION(ATTENTION) The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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⚠ WARNING(AVERTISSEMENT) Use of this ME Equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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⚠ WARNING(AVERTISSEMENT) Use of accessories, transducers and cables other than those specified or provided by DENTIS of this ME Equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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⚠ WARNING(AVERTISSEMENT) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm(12 inches) to any part of the S200, including cables specified by DENTIS. Otherwise, degradation of the performance of this equipment could result.

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⚠ WARNING(AVERTISSEMENT) The instructions given in this document must be followed when handling the product. Failure to do so may endanger the safety of the installers or users.

As well as specific information on operating the entire product and conducting preventive maintenance, are provided in the USER'S MANUAL. For further information, please contact our sales network or our local network.

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⚠ WARNING(AVERTISSEMENT) The electrical connections must be performed by a qualified technician only.

The electrical installation must be planned, performed and inspected by electrical engineers.

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No	Location
①	Marking on the outside of ME EQUIPMENT
②	Marking on the inside of ME EQUIPMENT
③	Marking on the controls of ME EQUIPMENT
④	Marking on the packing label of ME EQUIPMENT
⑤	Marking on the manual of ME EQUIPMENT
⑥	Marking on the label of ME EQUIPMENT

4. Introduction

4.1 Intended use

- Minor Surgical Luminaire(S200) is intended to be used to provide visible illumination of the surgical area or the patient during surgical operations, diagnosis and treatment.
- Classification under the provision of 93/42/EEC(MDR) : Class I
- Minor Surgical Luminaire(S200) is classified as a Class I device.
- Classification under the provision of FDA (U.S. Food and Drug Administration) : Class I
- Minor Surgical Luminaire(S200) is classified as a Class I device.

- Form of protection against electric shock : Class I
 - Minor Surgical Luminaire(S200) is classified as Class I
- Degree of protection against flammability
 - Minor Surgical Luminaire(S200) is classified as a device not suitable to be used in a potentially flammable environment.
 - Do not use near flammable materials.
- Method(s) of sterilization or disinfection recommended by the manufacturer.
 - The SIDE HANDLE should be cleaned with cloths, and the MAIN HANDLE should be sterilized with sterilizer regularly to prevent disinfection.
- Mode of operation
 - Classification of Minor Surgical Luminaire(S200) : continuous operation.

4.2 General description

- The user must ensure that that the device works properly and is in a satisfactory condition before each use.
- This DENTIS device is intended only for use in the field of medical. It is impermissible to use the device for a purpose for which it was not intended.
- "Proper Use" includes following all the instructions for use and ensuring that all inspections and service tasks are performed.
- Apply and meet the overarching guidelines and/or national laws, national regulations and the rules of technology for medical devices applicable for startup and use of the device for the intended purpose.
- The user must observe the following :
 - Only use properly operating equipment.
 - Protect himself or herself and third parties from danger.
 - Avoid contamination from the device.
- During use, the following national regulations must be observed:
 - The applicable health and safety regulations.
 - The applicable accident prevention regulations.
- To ensure that device maintains their value and are always ready for use, they must be serviced once a year as recommended.
- Before using the products, You must receive training by authorized person of the DENTIS.
- The safety checks must be performed every year.
- Repair and service of the device is authorized only to those who meet the requirements below:
 - Technicians of authorized dealers specially trained by DENTIS.
 - The trained technicians of DENTIS branches.

4.3 Environmental requirement

- | | |
|--|---|
| <ul style="list-style-type: none"> · Conditions of the usage environment <ul style="list-style-type: none"> - Temperature : 0 - 40°C - Relative Humidity : 30 - 90 %RH - Atmospheric pressure: 0-2,000 m (106-80 Kpa) | <ul style="list-style-type: none"> · Conditions of storage and transportation environment <ul style="list-style-type: none"> - Temperature : 0 - 40 °C - Relative Humidity : 0 - 80 % - Atmospheric pressure: 0 - 2,000 m (106 - 80 Kpa) |
|--|---|

4.4 Safety information

- Minor Surgical Luminaire(S200) as a medical device complies with the safety regulation EN/IEC 60601-1.
- Accessory equipment connected to the analog and digital interfaces must be certified according to the respective EN/IEC standards (e.g. EN/IEC 60601-2-41 for Particular requirements, EN/IEC 60601-1 for General requirements).

- Furthermore all configurations shall comply with the ME Equipment of Collateral Standard EN/IEC 60601-1-2. Everybody who connects additional equipment to the signal input part or signal output part configures a ME Equipment, and is therefore responsible that the ME Equipment complies with the requirements of standard EN/IEC 60601-1.

- If in doubt, consult the technical service department or your local representative.

- For EU Countries.

Europe Representative
KTR Europe GmbH
Merqenthalerallee 77, 65760 Eschborn, Germany

- Safety is everyone's obligation and responsibility.
- The safe use of this product is related to everyone such as installer, user, operator and equipment's manager.
- It must read and learn this user's manual is compulsory before installation, using, cleaning, fixing or operation of this product or its accessories. Pay particular attention and be familiar with warning symbols about safety.
- If do not follow safety direction of this manual, you can get injured or accident when you operate this product. After read carefully and understand this manual, use this product.
- This manual is in keep a place where you can find easily.

4.5 Warranty regulation

- DENTIS warrants all products against defects in materials or workmanship for two year from time of delivery.
- DENTIS's sole obligation under the warranty is to provide parts for the repair, or at its option, to provide the replacement product (excluding labor).
- The buyer shall have no other remedy. All special, incidental, and coincidental damages are excluded.
- Written notice of breach must be given to DENTIS within the warranty period.
- The warranty does not cover damage resulting from improper installation or maintenance, accident or misuse.
- The warranty does not cover damage resulting from the use of cleaning, disinfecting or sterilization chemicals and processes.
- The Failure to follow instructions provided in the DENTIS Instructions for Use (operation and maintenance instructions) may void the warranty.
- LED PCB Ass'y is covered under 60,000 hrs warranty.

4.6 Electromagnetic Compatibility

4.6.1 Emission

This ME equipment is intended for use in Professional healthcare facility environment.		
Emission test	Compliance	Guidance
Conducted Disturbance CISPR 11(EN 55011)	Complies (Group 1, Class A)	The S200 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated Disturbance CISPR 11(EN 55011)		
Harmonic current IEC 61000-3-2	Complies	The S200 is suitable for use in all establishments other than domestic premises and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations and flicker IEC 61000-3-3	Complies	

4.6.2 Immunity

This ME equipment is intended for use in Professional healthcare facility environment.					
Immunity test	EN 60601-1-2:2015		Compliance		
Electrostatic Discharge(ESD) IEC 61000-4-2	Direct: ± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air Indirect: ± 8 kV HCP/VCP		Complies		
Radio Frequency Electromagnetic Fields IEC 61000-4-3	3 V/m @ 80 MHz ~ 2.7 GHz 80 % AM at 1 kHz		Complies		
Proximity fields from RF wireless communications equipment IEC 61000-4-3	Frequency (MHz)	Modulation	Immunity Level (V/m)		
	385	**Pulse Modulation: 18 Hz	27		
	450	*FM ± 5Hz deviation: 1 kHz sine	28		
	710 745 780	**Pulse Modulation: 217 Hz	9		
	810 870 930	**Pulse Modulation: 18 Hz	28		
	1720 1845 1970	**Pulse Modulation: 217 Hz	28		
	2450	**Pulse Modulation: 217 Hz	28		
	5240 5500 5785	**Pulse Modulation: 217 Hz	9		
	** The carrier shall be modulated using a 50 % duty cycle square wave signal. * As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.			Complies	
	Fast Transients IEC 61000-4-4	Voltage	AC/DC power ports		Signal ports
Test voltage		± 2 kV	± 1 kV		
- 100 kHz repetition frequency					
Surges IEC 61000-4-5	Voltage	Power lines			
	Test voltage	Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV			
RF Continuous Conducted IEC 61000-4-6	3 V @ 0.15 MHz ~ 80 MHz 6 V @ in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		Complies		
Power Frequency Magnetic Fields IEC 61000-4-8	30 A/m @ 50 Hz or 60 Hz		Complies		
Voltage Dips, Interruptions, and Variations IEC 61000-4-11	<ul style="list-style-type: none"> • Voltage Dips 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° • Voltage Interruptions 0 % UT; 250/300 cycle • Voltage Variations 				Complies
	Frequency (Hz)	Ranges			
	50	49, 50, 51			
	60	59, 60, 61			

5. LIGHTHEAD specification

5.1 Product colors : BLUE, IVORY, WHITE

5.2 Technical specification (In accordance with EN/IEC 60601-2-41)

· STANDARD SPECIFICATION

Specifications	Unit	Luvis S200		
		60,000 lx	80,000 lx	100,000 lx
Normal Mode		60,000 lx	80,000 lx	100,000 lx
Central illuminance Ec (@1m)	lx	60,000	80,000	100,000
Electronic adjustment range	%	8.3~100	6.3~100	5~100
Light field diameter (d10)	cm	18	18	18
Light field diameter (d50)	cm	9.5	9.5	9.5
Light field uniformity (d50/d10)	N/A	> 0.5	> 0.5	> 0.5
Depth of Illumination (60%)	cm	78	78	78
Color temperature (3 levels)	K	3,800 / 4,300 / 4,800		
Color rendering index (Ra)	N/A	90	90	90
Special color rendering index (R9)	N/A	-	-	-
Shadow dilution	With tube	%	100	100
	With one mask	%	-	-
	With two mask	%	48	48
	With one mask, at base of tube	%	-	-
	With two mask, at base of tube	%	46	46
Radiation energy	mW/m ² ·lx	3.26	3.32	3.32
Illumination energy(Ee)	W/m ²	204	265	319

* Optical values are measured with a tolerance of ±10%

5.3 Electrical specification (In accordance with EN/IEC 60601-1)

· SMPS TECHNICAL DATA

Content	Description
Input voltage	AC 100-240 V, 50/60 Hz
Output voltage	DC 33 V, 0.5 A

· RATING

Content	Rated Input	Remark
60,000 lx	AC 100-240 V, 50/60 Hz, 27-43 VA	-
80,000 lx	AC 100-240 V, 50/60 Hz, 34-52 VA	-
100,000 lx	AC 100-240 V, 50/60 Hz, 45-67 VA	-

6. USE

6.1 Product type

Dual connection type

Wall type

Mobile type

Single ceiling type

Dual ceiling type

5.4 Mechanical specification


Specifications	Length(mm)	Weight(kg)	Remark
LIGHT HEAD	346 × 348 × 98	3.5	+HEAD ARM
TENSION SPRING ARM	Φ32 × 650	2.50	-
WALL SECOND ARM	Φ34 × 707	2.04	-
DUAL CEILING SECOND ARM	Φ34 × 707	2.04	-
CEILING SECOND ARM	Φ32 × 695	1.61	-
CHAIR SECOND ARM	Φ32 × 777	4.15	-
MOBILE VERTICAL ARM	Φ31.8 × 1,693	3.05	+MOBILE VERTICAL ARM1 +MOBILE VERTICAL ARM2
MOBILE BASE	575 × 575	17.05	+CASTER 4EA + BASE WEIGHT 2EA
DUAL CONNECTOR	132 × 71	1.26	+ARM BOLT
WALL BRACKET	100 × 270	2.60	-
CEILING VERTICAL ARM	Φ42.7 × 733	4.82	-
DUAL CEILING VERTICAL ARM	Φ60.5 × 865	17.10	+DUAL BRACKET
CEILING COVER & BRACKET	Φ500 × 50	1.15	-
CEILING MOUNT	-	13.36	GUIDE BRACKET + FIX BRACKET + STUD BOLT

5.5 Other characteristic


Specifications	Luvis S200	Remark
Protection against electric shock	Class I Protection	-
Protection against harmful ingress of water or particulate matter	Ordinary	Head part: IP43
Method of sterilization	See the USER'S MANUAL 8. Sterilization and Cleaning	-

6.2 Product description

· LIGHTHEAD


· Pattern Size


UNIT : mm

· Controller Panel


· Depth

6.3 Performance characteristic

· MOBILE

(unit : mm)


⚠ While moving the mobile type device for use, please do not step on the power adapter and do not place it on the ground.

⚠ While moving the mobile type device for use, please refer to the instruction and the picture on the step 10 of 9. OPERATION.

⚠ Do not push the mobile type device when wheels are locked.


· DUAL CONNECTION

(unit : mm)


· WALL

(unit : mm)


· SINGLE CEILING

(unit : mm)

· DUAL CEILING

(unit : mm)

⚠ If this device is operated beyond the indicated angles, it can be damaged.

⚠ This image was constructed to demonstrate the product's operating range. The image may differ from the actual product.


6.4 Detachable part

MAIN HANDLE	2ea	
-------------	-----	--

⚠ Before using MAIN HANDLE, please conduct a sterilization treatment.


7. Labelling requirement

7.1 MARKING PLATE

NO.	LABEL NAME	LABEL IMAGE	SIZE(mm)				
①	PRODUCT LABEL		84 X 54				
②	No pushing		20 X 20				
③	Do not use this device in a bath-tub, shower or water-filled reservoir		20 X 20				
④	MAIN HANDLE LABEL (On the plastic cover of MAIN HANDLE)	<table border="1" data-bbox="706 1619 1080 1780"> <tr> <td>PRODUCT : MAIN HANDLE</td> </tr> <tr> <td>LOT NO:</td> </tr> <tr> <td>Manufacturer : DENTIS CO., LTD.</td> </tr> <tr> <td>99,Seongseoseo-ro, Dalseo-gu, Daegu, Korea</td> </tr> </table>	PRODUCT : MAIN HANDLE	LOT NO:	Manufacturer : DENTIS CO., LTD.	99,Seongseoseo-ro, Dalseo-gu, Daegu, Korea	100 X 45
PRODUCT : MAIN HANDLE							
LOT NO:							
Manufacturer : DENTIS CO., LTD.							
99,Seongseoseo-ro, Dalseo-gu, Daegu, Korea							

7.2 CAUTION AND INFORMATION MARK

⚠ Below is the reference for label positions. Please refer to 7.1 Marking Plate for matching the indicated numbers to their name plates.

8. Sterilization and Cleaning

⚠ Cleaning and sterilization procedures vary from facility to facility and therefore it is not possible to specify a single procedure. The user should contact a specialist in the hospital and apply the recommended procedures according to the product. If you have any questions about the compatibility of your medication, please contact the DENTIS Customer Center.

The treatment level required for cleaning and disinfection of the MAIN HANDLE is a low level of disinfection. It is classified as non-critical devices with a low risk of infection, except for the use of sterilizable handles.

Contact your place of purchase in case of failure or damage.

⚠ The use of cleaning agents containing the following substances is not permitted.
- High concentrations of organic and inorganic acids and chlorinated hydrocarbons

8.1 Cleaning the equipment

⚠ Before cleaning turn off the power and wait until LIGHTHEAD is sufficiently cooled.
Wear gloves when cleaning and disinfecting.
Do not spray cleaning and disinfectant directly into the LIGHTHEAD.

- Make sure that the system is powered off and the LIGHTHEAD is sufficiently cooled.
- Prepare cleaning agents and disinfectants. For cleaning agents, follow the manufacturer's instructions for use.
- Moisten a clean, lint-free, soft cloth with cleaning solution and remove excess moisture.
- Clean the area using a cloth soaked with cleaning agent.
- Clean the device using a cloth soaked with clean water.
- Wipe all surfaces with a clean, dry cloth to remove any residue.
- Make sure all cleaning and disinfectant residues have been removed before using the light.

8.2 MAIN HANDLE AUTOMATIC CLEANING

- After the treatment, remove the MAIN HANDLE from the dental lamp.
- Remove visible foreign substances using a soft brush and cloth in warm water flowing through the MAIN HANDLE for at least 1 minute.
- Place the MAIN HANDLE in the cleaning equipment and carry out the cleaning using an enzyme detergent according to the following cleaning conditions.

No.	Condition	Temperature(°C)	Minimum time(Min)
1	Pre-cleaning	20-40	1
2	Cleaning	20-40	5
3	Cleaning	20-40	2
4	Rinsing	20-40	10
5	Rinsing	20-40	10
6	Dry	50-70	90

- Remove the MAIN HANDLE from the cleaning equipment and ensure that the cleaning is complete.
- Check for foreign substances left inside and outside the MAIN HANDLE, if it is necessary, repeat the cleaning process
- Protect the cleaned and dried MAIN HANDLE from re-contamination.

8.3 MAIN HANDLE MANUAL CLEANING

- After the treatment, remove the MAIN HANDLE from the dental lamp.
- Immerse the MAIN HANDLE in the washing tank for at least 5 minutes in the diluted enzyme detergent. (Follow the manufacturer's instructions for use of the cleaner.)
- To prevent blood or foreign substances from drying out, soak the separated MAIN HANDLE in distilled water or tap water for 15 minutes to remove any foreign substances, then clean it with a soft brush and a lint-free cloth.
- Check inside and outside of the MAIN HANDLE for debris, and repeat cleaning if any debris remains.
- Rinse thoroughly with clean water and wipe clean with a lint-free cloth to dry.
- Protect the cleaned and dried main handle from re-contamination




8.4 MAIN HANDLE Sterilization




⚠ Do not allow foreign substances to enter the MAIN HANDLE during sterilization.
MAIN HANDLE is guaranteed for up to 100 sterilization cycles when the same sterilization conditions are applied.
In the case of MAIN HANDLE, sterilization may cause wear, cracks, and discoloration. If these signs are detected, stop using the handle and replace it.
The MAIN HANDLE being stored must be kept clean and sterile before use.
If you are not wearing a sterile suit, be careful and do not touch the sterilized MAIN HANDLE.

9. Operation

- Sterilize the cleaned MAIN HANDLE according to the process below.
- Make sure the MAIN HANDLE is clean, wrap the HANDLE in sterile packaging (double packaging or equivalent) and seal it.
- When putting the MAIN HANDLE into the sterilizer, make sure that the hole is facing down to allow water to flow down.
- Place the MAIN HANDLE in the Pressurized Steam Sterilizer (AUTO CLAVE) and proceed with the sterilization cycle according to the sterilization conditions. Using chemical sterilization is prohibited
- Follow the sterilizer manufacturer's instructions and sterilization conditions of your country.
- When storing after sterilization, store in a tray in a sterile packaging material.(room temperature storage)

Sterilization Condition	
<ul style="list-style-type: none"> · Operating Condition <ul style="list-style-type: none"> - Temperature : 132 °C - Pressure : 160 ± 50kPa (1.6 ±0.5 kgf/cm²) - Time : 10 min - Dry : 16 min 	 <p>Autoclave</p>
Sterilization Condition	
<ul style="list-style-type: none"> · Sterile Clothing <ul style="list-style-type: none"> - KIMGUARD sterile clothing - Model : KC500 	

Step 1	Before using the device, check the maneuverability of BOTTOM ASSY S2 by moving it up and down, left and down.	
Step 2	As shown in the image next, please plug the power cord to the installed outlet. (DUAL CONNECTION, WALL, MOBILE) Plug the power cord to an outlet, shown and described in the picture. For a Mobile type, turn on the power switch, located at the bottom of the VERTICAL ARM's handle.	
	For a Single Ceiling / Dual Ceiling type, a wall switch, as shown in the picture, must be installed first. To use the device, turn on the switch.	
Step 3	The LED turns on when the power is applied.	
Step 4	Before and after using BOTTOM ASSY S2 HANDLE, please sterilize it with a specified disinfectant.	


Step 5	<p>Select the desired Color temperature (ORANGE, PINK, WHITE)</p> <ul style="list-style-type: none"> - Orange : 3,800K - Pink : 4,300K - White : 4,800K <p>Select the desired luminosity(5steps)</p>	
Step 6	Move the head part to the desired position.	
Step 7	<p>When MAIN HANDLE is rotated clockwise, Pattern size is 180mm</p> <p>When MAIN HANDLE is rotated counterclockwise, Pattern size is 300mm</p>	
Step 8	After using of device, turn off the STAND-BY button.	

	After using the device, MAIN HANDLE should be detached and sterilized.	
	⚠ Please pull out MAIN HANDLE by pressing the fixing key, marked with a red arrow.	
Step 9	⚠ Always use the sterilized MAIN HANDLE.	
	<p>⚠ Sterilization (Autoclave)- Service Condition</p> <p>Temperature: 132 °C</p> <p>Pressure: 160 ± 50kPa (1.6 ±0.5 kgf/cm²)</p> <p>Time: 10 min</p> <p>Dry: 16 min</p>	
Step 10	<p>The picture 1 demonstrates the method of moving the mobile type</p> <p>⚠ Check the wheel LOCK is OFF</p> <p>Storage method of device:</p> <p>Devices other than the mobile type must be positioned so that the user and patient can avoid any collisions with the devices.</p> <p>(Refer to the picture 1 and picture 2)</p>	

10. Adjustment and Maintenance

10.1 Adjustment

⚠ A yearly checkup on the intensity and Color temperature of the device at an authorized dealer is recommended.

If the head does not stay fixed and slides left or right,
Use a M5 wrench to adjust the tension of head joint.

- Clockwise : Higher tension
- Counterclockwise : Lower tension

⚠ If the joint is wrenched clockwise too tightly, the joint arm may get damaged.

⚠ Do not attempt to repair any damage or wear and tear on the device.
Please contact an authorized dealer for help.

10.2 Maintenance

Trouble 1	If the LED light on the head does not operate.	<ul style="list-style-type: none"> - Check if the power code is plugged in the outlet properly. - Check if the STAND-BY INDICATOR LED in front of the control panel is on and green. - Press the Control panel stand-by button and check that the LED turns on. (If green light does not light on, do not attempt to repair it yourself, contact your local authorized dealer.)
Trouble 2	If head is not fixed after handling, moves left and right.	 <ul style="list-style-type: none"> - Check that the verticality and horizontality of the installed device using circular leveler. As shown in the picture, use the M3 wrench to tighten it. If the condition does not improve, do not attempt to repair. Please contact your local authorized dealer for help.
Trouble 3	If there is an error in the pattern formation.	<p>Hold MAIN HANDLE and check if rotation is prospered. (Refer operation step) If the condition does not improve, do not attempt to repair. Please contact your local authorized dealer for help.</p>
Trouble 4	If the light shield cover is polluted.	<p>Cleanse it with designated chemicals (Alcohol, Ethanol). If the essential performances of the device (Intensity of Illumination, Color Temperature) are seriously impaired, please contact your local authorized dealer for repair.</p>

11. Disposal

- For environment and safety of human, wastes must be recycled or separated.
- The materials should be carefully separated.
- The electrical boards should be submitted to an appropriate recycling proceeding
- The cardboard box may be recycled with other paper products.


- Please contact the closest DENTIS branch or your supplier, if you have any questions about recycling of the device

12. Model designation

Model	ARM Axis	Mount Type	Colour	Intensity	Power supply
S200	2=2 Axis	0=Head Only 1=Ceiling Type 2=Dual Ceiling Type 3=Wall Type 4=Dual Connection Type 5=Mobile Type	0=Blue 1= Ivory 2=White	0=100,000	0 = PLUG 1=External SWITCH 2=POWER SWITCH


13. Manufacturer and Factory

Manufacturer, Factory : DENTIS MEDICAL DIVISION
6, Yuram-ro, Dong-gu, Daegu, Republic of Korea

14. Europe Representative

Location of KTR Europe GmbH. Mergenthalerallee 77, 65760 Eschborn, Germany

DENTIS

DENTIS Co., Ltd. (HEADQUARTERS) 87, Yuram-ro Dong-gu, Daegu, Republic of KOREA Tel +82-53-583-2804 Fax +82-53-583-2806
DENTIS MEDICAL DIVISION(FACTORY) 6, Yuram-ro, Dong-gu, Daegu, Korea Tel. +82-70-7730-4146 Email. luvis@luvis.co.kr

EL2SAP-PS080