Professional LED Light System  
for Medical Surgery

# LUVIS S300 USER'S MANUAL

---

**Minor Surgical Luminaire**

**DENTIS CO., LTD.**

**87, Yuram-ro, Dong-gu, Daegu, Republic of Korea**

Tel. +82-53-583-2804, FAX. +82-53-583-2806 [www.dentis.co.kr](http://www.dentis.co.kr) / [www.luviz.co.kr](http://www.luviz.co.kr)

## CONTENTS

1. Standards .....	5
2. Cautions and Warnings .....	6
3. Symbol(Symbole).....	14
4. Introduction .....	16
4.1 Intended use .....	16
4.2 General description .....	17
4.3 Environmental requirement.....	18
4.4 Safety information.....	18
4.5 Warranty regulation.....	19
4.6 Electromagnetic Compatibility .....	20
5. LIGHTHEAD specification.....	22
5.1 Technical specification (In accordance with EN/IEC 60601-2-41, EN/ISO 9680) .....	22
5.2 Electrical specification (In accordance with EN/IEC 60601-1) .....	24
5.3 Mechanical specification .....	24
5.4 Other characteristic.....	24
6. USE.....	26
6.1 Controlling the LIGHTHEAD with the ARM CONTROLLER .....	26
6.2 Light functions operated at the HANDLE CONTROLLER .....	28
6.3 Operation of the surgical light head using the REMOTE CONTROLLER .....	28
6.4 Positioning.....	29
6.5 Light field diameter .....	30
6.6 MAIN HANDLE .....	31

6.7 Product Type .....	32
6.8 Mobile Type.....	35
6.9 BATTERY PACK.....	37
7. Android application(Optional).....	38
7.1 Controlling the LIGHTHEAD with the Android application(optional).....	38
7.2 Download the LUVIS CONTROL APP.....	42
7.3. Installing and setting the LUVIS CONTROL APP.....	46
8. Performance characteristic.....	53
8.1 Single Ceiling Type .....	53
8.2 Dual Ceiling Type.....	54
8.3 Wall Type.....	55
8.4 Dual connection Type .....	56
8.5 Mobile Type.....	57
9. Cleaning and Sterilization .....	59
9.1 Cleaning the equipment.....	59
9.2 MAIN HANDLE AUTOMATIC CLEANING.....	60
9.3 MAIN HANDLE MANUAL CLEANING.....	60
9.4 MAIN HANDLE Sterilization .....	61
10. Maintenance.....	62
11. Troubleshooting.....	63
11.1 LIGHTHEAD.....	63
11.2 Android application(Optional) .....	63
12. Disposal.....	65

13. Model designation ..... 66

14. List of component ..... 67

15. Manufacturer and Factory ..... 70

16. Europe Representative ..... 71

**“WARNING: Modification of this equipment is not allowed”**

## 1. Standards

- Certification of DENTIS
  - EN ISO 13485:2016
  - Relevant EC Regulation: REGULATION (EU) 2017/745
- Applied Standards:
  - EN ISO 15223-1:2016, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
  - EN 1041:2008/A1:2013, Information supplied by the manufacturer with medical devices
  - EN ISO 13485:2016, Medical devices – Quality management systems – Requirements for regulatory purpose
  - EN ISO 14971:2019, Medical devices – Application of risk management to medical devices
  - EN 60601-1:2006+A2:2021, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
  - EN 60601-1-2:2015, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances – Requirements and tests
  - EN 60601-1-6:2010, Medical electrical equipment-Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability
  - EN ISO 7010:2019, Graphical symbols – Safety colors and safety signs-Registered safety signs
  - EN 62471:2008, Photo biological safety of lamps and lamps systems
  - EN 62366-1:2015, Medical devices – Application of usability engineering to medical devices
  - IEC 60601-2-41:2009/A1:2015, Medical electrical equipment – Part 2-41: Particular requirements for the basic safety and essential performance of surgical luminaires and luminaires for diagnosis
  - EN ISO 9680:2021, Dentistry – Operating lights
  - IEC 62133:2012, Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications
  - IEC 62304:2006/A1:2015, Medical device software, Software life-cycle processes

## 2. Cautions and Warnings

### CAUTION(ATTENTION)

This ME Equipment is intended only for use in the professional healthcare facility environment, e.g. public and private hospitals, specialized medical offices, etc. This ME Equipment is intended for direct use on operation theatres and nearby HF surgical equipments, where the intensity of EM DISTURBANCES complies with the applicable standards. For further instructions please follow the chapter “Electromagnetic compatibility” of the User's Manual.

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### CAUTION(ATTENTION)

If this ME Equipment is lost or degraded the performance due to EM DISTURBANCES, result in improper operation and degradation of the performance.

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### CAUTION(ATTENTION)

All cables and maximum lengths of coaxial cables that are replaceable by the DENTIS and that are likely to affect compliance of the this ME Equipment with the requirements of EMC(Electro-Magnetic Compatibility).

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### CAUTION(ATTENTION)

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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### WARNING(AVERTISSEMENT)

Use of this ME Equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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### WARNING(AVERTISSEMENT)

Use of accessories, transducers and cables other than those specified or provided by DENTIS of this ME Equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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### WARNING(AVERTISSEMENT)

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm(12 inches) to any part of the S300, including cables specified by DENTIS. Otherwise, degradation of the performance of this equipment could result.

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### **WARNING(AVERTISSEMENT)**

The instructions given in this document must be followed when handling the product. Failure to do so may endanger the safety of the installers or users.

As well as specific information on operating the entire product and conducting preventive maintenance, are provided in the USER'S MANUAL. For further information, please contact our sales network or our local network.

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### **WARNING(AVERTISSEMENT)**

The electrical connections must be performed by a qualified technician only.

The Lithium Battery replacement must be performed by a qualified technician only. The electrical installation must be planned, performed and inspected by electrical engineers.

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### **WARNING(AVERTISSEMENT)**

The LIGHTHEAD is designed to operate using a AC 100–240V 50/60Hz. Higher or lower voltages may affect the light intensity and operating life of the LEDs.

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### **WARNING(AVERTISSEMENT)**

Damaged wire insulation may result in the risk of electric shock.

To protect against the risk of electrocution, connect the power cables carefully.

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### **WARNING(AVERTISSEMENT)**

Major Surgical Luminaire(LUVIS S300) is Class I equipment. In order to avoid the risk of an electric shock, the equipment must be connected to a mains supply with PE(Protective Earth).

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**WARNING(AVERTISSEMENT)**

Take care when handling the circuit boards: these boards are supplied in an electrostatic envelope and must be handled with great care.

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**WARNING(AVERTISSEMENT)**

Do not look directly into light source(LED).

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**WARNING(AVERTISSEMENT)**

The LIGHTHEAD brakes are adjusted during installation. Like all mechanical parts, the brakes are subject to wear.

Read just the brakes if the LIGHTHEAD no longer remains steady in any position.

Check the condition of the mounting surface.

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**WARNING(AVERTISSEMENT)**

The operation and safety of the device may be affected by the removal of certain components during servicing operations.

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**WARNING(AVERTISSEMENT)**

All the information in this manual has been checked out carefully and discerned as accurate one at the time of publication.

However, DENTIS takes no responsibilities of the results caused by default, omission, or misuse of it.

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**WARNING(AVERTISSEMENT)**

DENTIS has rights to modify the product itself or specifications of the product without any prior notice, as well as rights not to renew that modification on this manual.

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### **WARNING(AVERTISSEMENT)**

Do not install any other software on your TABLET PC. Maintenance by software is not our responsibility. If there is an error in installing the APP, please back up the data and factory reset the TABLET PC and reinstall the APP.

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### **WARNING(AVERTISSEMENT)**

The LUVIS CONTROL APP is subject to change and upgrade without notice.

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### **WARNING(AVERTISSEMENT)**

To upgrade the version of the LUVIS CONTROL APP, please check the homepage for the latest version.

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### **WARNING(AVERTISSEMENT)**

Please notify the manufacturer of bugs and errors in LUVIS CONTROL APP.

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### **WARNING(AVERTISSEMENT)**

Customer service is not provided regarding the LUVIS CONTROL APP.

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### 3. Symbol(Symbole)

Symbol (Symbole)	Meaning (Signification)	Remark (Remarque)
	<p><b>CE Mark</b> The device bears the CE mark and complies with the requirements of Regulation (EU) MDR 2017/745 for medical device. - L'appareil porte le marquage CE et est conforme aux exigences du Règlement (UE) MDR 2017/745 relatif aux dispositifs médicaux.</p>	④, ⑤, ⑥
	<p>MEDICAL – GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH [ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012, CAN/CSA C22.2 No. 60601-1:14] - MATÉRIEL MÉDICAL – GÉNÉRALITÉS CONCERNANT LES CHOCS ÉLECTRIQUES, LES INCENDIES ET LES DANGERS MÉCANIQUES SEULEMENT CONFORMÉMENT AUX [ANSI/AAMIES 60601-1:2005/(R)2012 et A1:2009/(R)2012/(R)2012, CAN/CSA C222 No. 60-114]</p>	⑥
	<p>Recommendation - Recommandation</p>	⑤
	<p>Protective earth (ground) - Terre de protection</p>	②
	<p>Alternating current - Courant alternative</p>	②
	<p>Stand-by - Mode veille</p>	③
	<p>Do not throw away with general household waste - Ne pas jeter avec les déchets ménagers ordinaires</p>	⑤
	<p>Caution - ATTENTION</p>	⑤
	<p>Warning - AVERTISSEMENT</p>	⑤
	<p>Follow instructions for use - Suivre les instructions d'utilisation</p>	⑥
	<p>General mandatory action - Action obligatoire générale</p>	①
	<p>Pushing prohibited - Pousser interdit</p>	①

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## 4. Introduction

### 4.1 Intended use

- Minor Surgical Luminaire(S300) is intended to be used to provide visible illumination of the surgical area or the patient during surgical operations, diagnosis and treatment. Minor Surgical Luminaire(S300) is intended to be used to provide visible illumination of the dental light area or the patient during dental surgery, diagnosis and treatment.
- Classification under the provision of 93/42/EEC(MDR) : Class I
  - Minor Surgical Luminaire(S300) is classified as a Class I device.
- Classification under the provision of FDA (U.S. Food and Drug Administration) : Class II
  - Minor Surgical Luminaire(S300) is classified as a Class II(exempt) device.
- Form of protection against electric shock : Class I
  - Minor Surgical Luminaire(S300) is classified as Class I
- Degree of protection against flammability
  - Minor Surgical Luminaire(S300) is classified as a device not suitable to be used in a potentially flammable environment.
  - Do not use near flammable materials.
- Method(s) of sterilization or disinfection recommended by the manufacturer
  - The MAIN HANDLE should be sterilized with sterilizer regularly to prevent disinfection.
- Mode of operation
  - Classification of Minor Surgical Luminaire(S300) : continuous operation.

## 4.2 General description

- The user must ensure that that the device works properly and is in a satisfactory condition before each use.
- This DENTIS device is intended only for use in the field of medical. It is impermissible to use the device for a purpose for which it was not intended.
- “Proper Use” includes following all the instructions for use and ensuring that all inspections and service tasks are performed.
- Apply and meet the overarching guidelines and/or national laws, national regulations and the rules of technology for medical devices applicable for startup and use of the device for the intended purpose.
- The user must observe the following:
  - Only use properly operating equipment.
  - Protect himself or herself and third parties from danger.
  - Avoid contamination from the device.
- During use, the following national regulations must be observed:
  - The applicable health and safety regulations.
  - The applicable accident prevention regulations.
- To ensure that device maintains their value and are always ready for use, they must be serviced once a year as recommended.
- Before using the products, You must receive training by authorized person of the DENTIS.
- The safety checks must be performed every year.
- Repair and service of the device is authorized only to those who meet the requirements below:
  - Technicians of authorized dealers specially trained by DENTIS.
  - The trained technicians of DENTIS branches.

### 4.3 Environmental requirement

- Conditions of the usage environment
- Temperature : 0 - 40 °C
- Relative Humidity : 30 - 90 %RH
- Atmospheric pressure: 0-2,000 m (106–80 Kpa)
- Conditions of storage and transportation environment
- Temperature : 0 - 40 °C
- Relative Humidity : 0 - 80 %RH
- Atmospheric pressure: 0-2,000 m (106–80 Kpa)

### 4.4 Safety information

- Minor Surgical Luminaire(S300) as a medical device complies with the safety regulation EN/IEC 60601-1.
- Accessory equipment connected to the analog and digital interfaces must be certified according to the respective EN/IEC standards (e.g. EN/IEC 60601-2-41 for Particular requirements, EN/IEC 60601-1 for General requirements).
- Furthermore all configurations shall comply with the ME Equipment of Collateral Standard EN/IEC 60601-1-2. Everybody who connects additional equipment to the signal input part or signal output part configures a ME Equipment, and is therefore responsible that the ME Equipment complies with the requirements of standard EN/IEC 60601-1.
- If in doubt, consult the technical service department or your local representative.
- For EU Countries
- The following mark, the name & address of the EU Representative shows compliance of the instrument with Directive Council Directive 93/42/EEC of 14 June 1993 as amended by Directive 2007/47/EC concerning medical devices.

- Europe Representative  
KTR Europe GmbH  
Mergenthalerallee 77, 65760  
Eschborn, Germany
- Safety is everyone's obligation and responsibility.
- The safe use of this product is related to everyone such as installer, user, operator and equipment's manager.
- It must read and learn this user's manual is compulsory before installation, using, cleaning, fixing or operation of this product or its accessories. Pay particular attention and be familiar with warning symbols about safety.
- If do not follow safety direction of this manual, you can get injured or accident when you operate this product. After read carefully and understand this manual, use this product.
- This manual is in keep a place where you can find easily.

#### 4.5 Warranty regulation

- DENTIS warrants all products against defects in materials or workmanship for two year from time of delivery.
- DENTIS's sole obligation under the warranty is to provide parts for the repair, or at its option, to provide the replacement product (excluding labor).
- The buyer shall have no other remedy. All special, incidental, and coincidental damages are excluded.
- Written notice of breach must be given to DENTIS within the warranty period.
- The warranty does not cover damage resulting from improper installation or maintenance, accident or misuse.
- The warranty does not cover damage resulting from the use of cleaning, disinfecting or sterilization chemicals and processes.
- The Failure to follow instructions provided in the DENTIS Instructions for Use (operation and maintenance instructions) may void the warranty.
- LED PCB Ass'y is covered under 60,000 hours warranty.

## 4.6 Electromagnetic Compatibility

### 4.6.1 Emissions

This ME equipment is intended for use in Professional healthcare facility environment.		
Emission test	Compliance	Guidance
Conducted Disturbance CISPR 11(EN 55011)	Complies (Group 1, Class A)	The ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated Disturbance CISPR 11(EN 55011)		
Harmonic current IEC 61000-3-2	Complies	The ME equipment is suitable for use in all establishments other than domestic premises and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations and flicker IEC 61000-3-3	Complies	

### 4.6.2 Immunity

This ME equipment is intended for use in Professional healthcare facility environment.				
Immunity test	EN 60601-1-2:2015		Compliance	
Electrostatic Discharge( ESD) IEC 61000-4-2	Direct: $\pm 8$ kV Contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV Air Indirect: $\pm 8$ kV HCP/VCP		Complies	
Radio Frequency Electromagnetic Fields IEC 61000-4-3	3 V/m @ 80 MHz ~ 2.7 GHz 80 % AM at 1 kHz		Complies	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	<b>Frequency (MHz)</b>	<b>Modulation</b>	<b>Immunity Level (V/m)</b>	
	385	**Pulse Modulation: 18 Hz	27	
	450	*FM $\pm 5$ Hz deviation: 1 kHz sine	28	
	710 745 780	**Pulse Modulation: 217 Hz	9	
	810 870 930	**Pulse Modulation: 18 Hz	28	
	1 720 1 845 1 970	**Pulse Modulation: 217 Hz	28	
	2 450	**Pulse Modulation: 217 Hz	28	
				Complies

	<table border="1"> <tr> <td>5 240</td> <td rowspan="3">**Pulse Modulation: 217 Hz</td> <td rowspan="3">9</td> </tr> <tr> <td>5 500</td> </tr> <tr> <td>5 785</td> </tr> </table> <p>** The carrier shall be modulated using a 50 % duty cycle square wave signal. * As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.</p>	5 240	**Pulse Modulation: 217 Hz	9	5 500	5 785		
5 240	**Pulse Modulation: 217 Hz	9						
5 500								
5 785								
Fast Transients IEC 61000-4-4	<table border="1"> <tr> <td>Voltage</td> <td>AC/DC power ports</td> <td>Signal ports</td> </tr> <tr> <td>Test voltage</td> <td>± 2 kV</td> <td>± 1 kV</td> </tr> </table> <p>- 100 kHz repetition frequency</p>	Voltage	AC/DC power ports	Signal ports	Test voltage	± 2 kV	± 1 kV	Complies
Voltage	AC/DC power ports	Signal ports						
Test voltage	± 2 kV	± 1 kV						
Surges IEC 61000-4-5	<table border="1"> <tr> <td>Voltage</td> <td>Power lines</td> </tr> <tr> <td>Test voltage</td> <td>Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV</td> </tr> </table>	Voltage	Power lines	Test voltage	Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV	Complies		
Voltage	Power lines							
Test voltage	Line to Line : ± 0.5 kV, ± 1 kV Line to ground: ± 0.5 kV, ± 1 kV, ± 2 kV							
RF Continuous Conducted IEC 61000-4-6	<p>3 V @ 0.15 MHz ~ 80 MHz 6 V @ in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz</p>	Complies						
Power Frequency Magnetic Fields IEC 61000-4-8	<p>30 A/m @ 50 Hz or 60 Hz</p>	Complies						
Voltage Dips, Interruptions, and Variations IEC 61000-4-11	<ul style="list-style-type: none"> <li>• Voltage Dips 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°</li> <li>• Voltage Interruptions 0 % UT; 250/300 cycle</li> <li>• Voltage Variations</li> </ul> <table border="1"> <tr> <td>Frequency (Hz)</td> <td>Ranges</td> </tr> <tr> <td>50</td> <td>49, 50, 51</td> </tr> <tr> <td>60</td> <td>59, 60, 61</td> </tr> </table>	Frequency (Hz)	Ranges	50	49, 50, 51	60	59, 60, 61	Complies
Frequency (Hz)	Ranges							
50	49, 50, 51							
60	59, 60, 61							

## 5. LIGHTHEAD specification

### 5.1 Technical specification (In accordance with EN/IEC 60601-2-41, EN/ISO 9680)

- STANDARD SPECIFICATION

Specifications		Unit	LUVIS S300	Remark
<b>Surgical Mode</b>				
Central illuminance $E_c$ (@1m)		lx	130,000	@4,800K
Electronic adjustment range		%	30~100	
Light field diameter (d10)		cm	15~25	
Light field diameter (d50)		cm	9~15	
Light field uniformity (d50/d10)		N/A	> 0.55	
Depth of Illumination (60%)		cm	75	
Color temperature (3 levels)		K	3,800/4,300/4,800	$\pm 300K$
Color rendering index (Ra)		N/A	96	$\pm 3$
Special color rendering index (R9)		N/A	95	$\pm 3$
Shadow dilution	With tube	%	100	
	With one mask	%	-	
	With two mask	%	50	
	With one mask, at base of tube	%	-	
	With two mask, at base of tube	%	50	
Radiation energy		$mW/m^2 \cdot lx$	3.62	
Illumination energy( $E_e$ )		$W/m^2$	454.5	

\* Optical values are measured with a tolerance of  $\pm 10\%$

- STANDARD SPECIFICATION(DENTAL MODE & EXAMINATION MODE)

Specifications	Unit	LUVIS S300	Remark
<b>Dental Mode</b>			
Central illuminance $E_c$ (@1m)	lx	5,000~30,000	
Electronic adjustment range	%	16~100 (3 level)	
Illuminance in patient's eyes	lx	<1,200	
Hard shadow	mm	12	
Focal Distance	cm	100	
Light field diameter (d10)	cm	15.0 X9.0	Oval shape
Light field diameter (d50)	cm	10.0 X6.0	Oval shape
Light field uniformity (d50/d10)	cm	> 0.65	
Color temperature	K	4,800	$\pm 300K$
Color rendering index (Ra)	N/A	95	$\pm 3$
<b>Examination Mode</b>			
Central illuminance $E_c$ (@1m)	lx	30,000~60,000	
Electronic adjustment range	%	50~100 (3 level)	
Light field diameter (d10)	cm	15	
Light field diameter (d50)	cm	9	
Light field uniformity (d50/d10)	cm	> 0.55	
Color temperature	K	3,800/4,300/4,800	$\pm 300 K$
Color rendering index (Ra)	N/A	96	$\pm 3$
Radiation energy	$mW/m^2 \cdot lx$	3.75	
Illumination energy( $E_e$ )	$W/m^2$	230.9	

\* Optical values are measured with a tolerance of  $\pm 10\%$

## 5.2 Electrical specification (In accordance with EN/IEC 60601-1)

- **SMPS TECHNICAL DATA**
- **MODEL : SMPS MODULE(Manufacturer : DENTIS CO., LTD.)**

Content	Single LIGHTHEAD	Dual LIGHTHEAD
Input	AC 100-240 V, 50/60 Hz	AC 100-240 V, 50/60 Hz x 2
Output	DC 24 V, 8.3A	DC 24 V, 8.3A x 2

- **RATING**

Content	Rating	Remark
S300 LIGHTHEAD	AC 100–240 V, 50/60 Hz, 84VA	with BATTERY PACK
	AC 100–240 V, 50/60 Hz, 61VA	without BATTERY PACK
BATTERY PACK	DC 14.8 V, 26 Ah	Lithium-ion Battery

## 5.3 Mechanical specification

Item	Length(mm)	Weight(kg)	Remark
LIGHT HEAD	422×409	4.9	+HEAD ARM
TENSION SPRING ARM	Φ32×650	2.50	-
WALL SECOND ARM	Φ34×707	2.04	-
DUAL CEILING SECOND ARM	Φ34×707	2.04	
CEILING SECOND ARM	Φ32×695	1.61	-
CHAIR SECOND ARM	Φ32×777	4.15	-
MOBILE VERTICAL ARM	Φ31.8×1693	3.05	+MOBILE VERTICAL ARM1+MOBILE VERTICAL ARM2
MOBILE BASE	600×586×197	20.5	+CASTER 4EA+BASE +SMPS CASE
DUAL CONNECTOR	132×71	1.26	+ARM BOLT

Professional LED Surgical Luminaire for All of Surgery Application

WALL BRACKET	100×270	2.60	-
CEILING VERTICAL ARM	Φ42.7×733	4.82	-
DUAL CEILING VERTICAL ARM	Φ60.5×865	17.10	+DUAL BRACKET
CEILING COVER & BRACKET	Φ500×50	1.15	-
CEILING MOUNT	-	13.36	GUIDE BRACKET + FIX BRACKET + STUD BOLT


#### 5.4 Other characteristic

Specifications	LUVIS S300	Remark
Protection against electric shock	Class I Protection	
Protection against harmful ingress of water or particulate matter	Ordinary	Head part: IP43
Method of sterilization	See the USER'S MANUAL	

## 6. USE

### 6.1 Controlling the LIGHTHEAD with the ARM CONTROLLER


#### 6.1.1 MAIN

#### FUNCTION

LUVIS S300		
	MAIN HEAD	<ul style="list-style-type: none"> <li>MAIN HEAD</li> </ul>
	SATELLITE HEAD	<ul style="list-style-type: none"> <li>SATELLITE HEAD</li> </ul>
	Synchronization	<ul style="list-style-type: none"> <li>Synchronizing LIGHTHEADs with each other.</li> <li>The setting is automatically applied to the Dual LIGHTHEAD.</li> </ul>
	Setting	<ul style="list-style-type: none"> <li>LIGHTHEAD setting</li> </ul>
	On/Off (Standby)	<ul style="list-style-type: none"> <li>On/Off button.</li> <li>The light turns on, at the last illumination level stored in its memory.</li> </ul>
	Level adjustment for the selected function	<ul style="list-style-type: none"> <li>Three light field diameter levels</li> <li>Three illumination levels</li> <li>Three color temperature levels 3,800K / 4,300K / 4,800K</li> </ul>
	Maximum mode	<ul style="list-style-type: none"> <li>Illumination : Level 3</li> <li>Light field diameter : Level 3</li> </ul>
	Examination mode	<ul style="list-style-type: none"> <li>Examination light mode.</li> </ul>
	Dental mode	<ul style="list-style-type: none"> <li>Dental light mode.</li> </ul>


6.1.2 SETTING

FUNCTION


LUVIS S300		
	Operating room No.	<ul style="list-style-type: none"> <li>Operating room number setting.</li> </ul>
	VERSION	<ul style="list-style-type: none"> <li>Firmware version</li> </ul>
	LIGHTHEAD	<ul style="list-style-type: none"> <li>Select MAIN HEAD or SATELLITE HEAD.</li> <li>If Single LIGHT, Select MAIN HEAD.</li> </ul>
	Number key pad	<ul style="list-style-type: none"> <li>Input the operating room number.</li> </ul>
	Page exiting	<ul style="list-style-type: none"> <li>Exiting from the setting page.</li> </ul>
	Backspace	<ul style="list-style-type: none"> <li>Erasing the operating room number.</li> </ul>
	Save	<ul style="list-style-type: none"> <li>Saving the setting value.</li> </ul>

### 6.2 Light functions operated at the HANDLE CONTROLLER

No.	Function	Description
①	Illumination level button	<ul style="list-style-type: none"> <li>The illumination level can be changed (Rotation) 1 → 2 → 3 → 1 → 2...</li> </ul>
②	Light field diameter level up button	<ul style="list-style-type: none"> <li>The light field diameter increases.</li> </ul>
③	Light field diameter level down button	<ul style="list-style-type: none"> <li>The light field diameter decreases.</li> </ul>

### 6.3 Operation of the surgical light head using the REMOTE CONTROLLER


**6.4 Positioning**

**Positioning the S300 LIGHTHEAD**


MAIN HANDLE

- Use the MAIN HANDLE to position the S300 LIGHTHEAD before the operation.
- This MAIN HANDLE can be removed for sterilization.

- Recommended distance : 1,000 mm

6.5 Light field diameter

**6.6 MAIN HANDLE**

**WARNING**

After each sterilization and before using the MAIN HANDLE again :

- Check MAIN HANDLE for cracks.
- Check that the MAIN HANDLE click into place correctly on the LIGHTHEAD.

**Attaching the MAIN HANDLE**


- Insert the MAIN HANDLE into the mount until it clicks into place

**Detaching the MAIN HANDLE**

- Press down on the PUSH BUTTON while removing the MAIN HANDLE.

**6.7 Product Type**


**<Mobile>**

**<Dual Ceiling>**

Professional LED Surgical Luminaire for All of Surgery Application

<Single Ceiling>

<Wall>

<Dual connection>


## 6.8 Mobile Type

### SMPS

#### WARNING

The mains plug should not be difficult to remove, place it in position (Means for isolation is disconnect mains plug)


- While moving the mobile type for use, please fold the spring arm as much as possible and hold the TENSION SPRING ARM by hands.
- And please make sure, that all CASTERS are moveable. (Brakes are loosened, no cables in the movement range of the stand)
- Storage method of device: The mobile type must be positioned so that the user and patient can avoid any collisions with the devices.

**SMPS and BATTERY PACK**

**WARNING**

The mains plug should not be difficult to remove, place it in position (Means for isolation is disconnect mains plug)

- While moving the mobile type for use, please fold the spring arm as much as possible and hold the TENSION SPRING ARM by hands.
- And please make sure, that all CASTERS are moveable. (Brakes are loosened, no cables in the movement range of the stand)
- Storage method of device: The mobile type must be positioned so that the user and patient can avoid any collisions with the devices.

**WARNING**

Do not lean and sit on this ME Equipment.  
Do not push when the wheels are fixed.

**WARNING**

Grab the SPRING ARM and VERTICAL ARM when moving this ME Equipment.

**WARNING**

Be careful of collision, obstruction and trip over the wire when moving this ME Equipment.

**WARNING**

When cleaning the floor using liquid, please move this ME Equipment elsewhere and keep it.

**6.9 BATTERY PACK**

**WARNING**

The lithium-ion BATTERY PACK is a consumable item and should be replaced every 12 months by service personnel.(It can vary depending on the number of use.)  
Check the BATTERY PACK capacity level by pushing the BATTERY CHECK BUTTON.

**WARNING**


The mains plug should not be difficult to remove, place it in position  
(Means for isolation is disconnect mains plug)

**WARNING**

After the first installation, you must fully charge the BATTERY PACK.

**The control panel of the BATTERY PACK CONTROLLER**

No.	Function	Description
①	Battery charging level indicator	<ul style="list-style-type: none"> <li>The charging level can be indicated 25%, 50%, 75%, 100%</li> </ul>
②	BATTERY CHECK BUTTON	<ul style="list-style-type: none"> <li>Press to check the charging level of the BATTERY PACK</li> </ul>

- The BATTERY PACK can be used for a minimum period of 3h with fully charged.(100%)

## 7. Android application(Optional)


### 7.1 Controlling the LIGHTHEAD with the Android application(optional)

#### Installing with the Android application

Refer to the LUVIS website for the installation of the application.

<http://www.luvis.co.kr/eng/>

#### 7.1.1 MAIN


#### FUNCTION

LUVIS S300		
	MAIN HEAD	<ul style="list-style-type: none"> <li>MAIN HEAD</li> </ul>
	SATELLITE HEAD	<ul style="list-style-type: none"> <li>SATELLITE HEAD</li> </ul>
	Synchronization	<ul style="list-style-type: none"> <li>Synchronizing LIGHTHEADs with each other.</li> <li>The setting is automatically applied to the Dual LIGHTHEAD.</li> </ul>
	Setting	<ul style="list-style-type: none"> <li>LIGHTHEAD setting</li> </ul>
	On/Off (Standby)	<ul style="list-style-type: none"> <li>On/Off button.</li> <li>The light turns on, at the last illumination level stored in its memory.</li> </ul>

	<p>Level adjustment for the selected function</p>	<ul style="list-style-type: none"> <li>• Three illumination levels</li> <li>• Three light field diameter levels</li> <li>• Three color temperature levels 3,800K / 4,300K / 4,800K</li> </ul>
	<p>Maximum mode</p>	<ul style="list-style-type: none"> <li>• Illumination : Level 3</li> <li>• Light field diameter : Level 3</li> </ul>
	<p>Examination mode</p>	<ul style="list-style-type: none"> <li>• Examination light mode.</li> </ul>
	<p>Dental mode</p>	<ul style="list-style-type: none"> <li>• Dental light mode.</li> </ul>

7.1.2 SETTING

FUNCTION

LUVIS S300		
	Operating room No.	<ul style="list-style-type: none"> <li>Operating room number setting.</li> </ul>
	LIGHTHEAD status indication	<ul style="list-style-type: none"> <li>Displaying the Wi-Fi communication status of the LIGHTHEAD and the TABLET.</li> </ul>
	VERSION	<ul style="list-style-type: none"> <li>Firmware version</li> </ul>
	Number key pad	<ul style="list-style-type: none"> <li>Input the operating room number.</li> </ul>
	Page Exiting	<ul style="list-style-type: none"> <li>Exiting from the setting page.</li> </ul>
	Backspace	<ul style="list-style-type: none"> <li>Erasing the operating room number.</li> </ul>
	Save	<ul style="list-style-type: none"> <li>Saving the setting value.</li> </ul>

**WARNING**

The LUVIS control application is recommended on Android 7.1. It may not work with other OS versions.

The LUVIS control application must be installed on 8-inch(1280 X 800 pixels) tablet PC. DENTIS is not responsible for any problems that arise when installing on an equipment other than 8-inch(1280x800 pixels) tablet PC.

**WARNING**

Do not install any other software on your tablet PC. Maintenance by software is not our responsibility. If there is an error in installing the app, please back up the data and factory reset the tablet PC and reinstall the app.

## 7.2 Download the LUVIS CONTROL APP

### 7.2.1 Download the LUVIS CONTROL APP using PC


- Connect the luvis website using PC and download the LUVIS CONTROL APP
- You need ID and password to download the LUVIS CONTROL APP
- If you are not registered as a member, please register as a member
- [http://luvis.co.kr/eng/content/04company/01\\_01.php](http://luvis.co.kr/eng/content/04company/01_01.php)


#### WARNING

If you can't connect to the website, contact your supplier or DENTIS.


- Press the button to download the application to a specific location.
- Connect TABLET and PC using the usb cable included with TABLET.

- Please put the downloaded LUVIS CONTROL APP on the TABLET.
- Download location : Tablet\Download


- Please check if LUVIS CONTROL APP downloading is completed in Installation files.
- File location : Tablet\Download


### 7.2.2 Download the LUVIS CONTROL APP using TABLET


- Connect the luvis website using TABLET and download the LUVIS CONTROL APP
- You need ID and password to download the LUVIS CONTROL APP
- If you are not registered as a member, please register as a member
- [http://luvis.co.kr/eng/content/04company/01\\_01.php](http://luvis.co.kr/eng/content/04company/01_01.php)

- Press the “APP Download” button to download the application to a specific location.
- Download location : Tablet\Download
- After downloading the LUVIS CONTROL APP, click the FORGET button to delete the connected Wi-Fi


- Please check if LUVIS CONTROL APP downloading is completed in “Installation files”.
- Download location : Tablet\Download

### 7.3. Installing and setting the LUVIS CONTROL APP

#### 7.3.1 Installing the LUVIS CONTROL APP

- Please go to the setting page of LIGHTHEAD Arm controller
- If you have existing LUVIS CONTROL APP installed, please delete it.
- Press Uninstall button
- Press OK button



- Follow the path below to see if there is an APP.



(My File → Installation files → LuvissApp(5INCH).apk)



- Make sure the version of the application is up to date.
- Press LuvissApp(5INCH).apk to install
- Press SETTINGS to enter the setting page

- You must allow “Allow from this source” in the setting  
(Setting → Security → Install unknown apps → My File)
- Activate “Allow from this source”

- Press the back button
- Press INSTALL
- Start the installation
- Press DONE


### 7.3.2 Setting the License Key of the LUVIS CONTROL APP


**WARNING**

If you lost license key, contact supplier or DENTIS


- Run the LUVIS CONTROL APP


- Please enter the dealer's email address and Password.


- Click certification button to generate Certification Key.

- Capture the generated certification key value and screen and send it to the DENTIS sales team using email.


- Enter the contents of License key in the email answered on the app.


- Click the OK button to go to the model selection step.

### 7.3.3 Initial setting of the LUVIS CONTROL APP


- Select the model
- Select the model and click the OK button
- Press the YES button

- Enter the same ROOM No. as LIGHTHEAD Room no.
- Press SAVE( Save ) button
- Press the YES button


- Press the OK button to connet LIGHTHEAD


### 7.3.4 Setting the Wi-Fi


- Slide from top to bottom to open the widget.
- Press and hold the Wi-Fi icon.
- Check the ROOM No. through the LIGHTHEAD's setting page.
- Select the same number of Wi-Fi as the ROOM No. in LIGHTHEAD.


- Enter Wi-Fi password(dentislight).
- Press the CONNECT.


- Check the Wi-Fi connection.

- Run the application.

- Press the On/Off button(  ) to operate the LIGHTHEAD.


- Check the Wi-Fi state


State	Description
	<ul style="list-style-type: none"> <li>• Before LIGHTHEAD selection</li> <li>• Disconnection Wi-Fi</li> </ul>
	<ul style="list-style-type: none"> <li>• Selection LIGHTHEAD</li> <li>• Disconnection Wi-Fi</li> </ul>
	<ul style="list-style-type: none"> <li>• Before LIGHTHEAD selection</li> <li>• Connection Wi-Fi</li> </ul>
	<ul style="list-style-type: none"> <li>• Selection LIGHTHEAD</li> <li>• Connection Wi-Fi</li> </ul>

## 8. Performance characteristic

### 8.1 Single Ceiling Type


<Unit : mm>


Professional LED Surgical Luminaire for All of Surgery Application

**8.2 Dual Ceiling Type**


(Unit : mm)


Professional LED Surgical Luminaire for All of Surgery Application


**8.3 Wall Type**


(Unit : mm)

**8.4 Dual connection Type**

<Unit : mm>

**WARNING**

If this device is operated beyond the indicated angles, it can be damaged. During positioning, eventual collisions between the LIGHTHEAD, SPRING ARMS and other devices must be avoided.

**WARNING**

This image was constructed to demonstrate the product's operating range. The image may differ from the actual product.

## 9. Cleaning and Sterilization

### Warning

Cleaning and sterilization procedures vary from facility to facility and therefore it is not possible to specify a single procedure. The user should contact a specialist in the hospital and apply the recommended procedures according to the product. If you have any questions about the compatibility of your medication, please contact the DENTIS Customer Center.

The treatment level required for cleaning and disinfection of the MAIN HANDLE is a low level of disinfection. It is classified as non-critical devices with a low risk of infection, except for the use of sterilizable handles.

Contact your place of purchase in case of failure or damage.

### Warning

The use of cleaning agents containing the following substances is not permitted.

- High concentrations of organic and inorganic acids and chlorinated hydrocarbons

### 9.1 Cleaning the equipment

### Warning

Before cleaning turn off the power and wait until LIGHTHEAD is sufficiently cooled.

Wear gloves when cleaning and disinfecting.

Do not spray cleaning and disinfectant directly into the LIGHTHEAD.

- Make sure that the system is powered off and the LIGHTHEAD is sufficiently cooled.
- Prepare cleaning agents and disinfectants. For cleaning agents, follow the manufacturer's instructions for use.
- Moisten a clean, lint-free, soft cloth with cleaning solution and remove excess moisture.
- Clean the area using a cloth soaked with cleaning agent.
- Clean the device using a cloth soaked with clean water.
- Wipe all surfaces with a clean, dry cloth to remove any residue.
- Make sure all cleaning and disinfectant residues have been removed before using the light.

## 9.2 MAIN HANDLE AUTOMATIC CLEANING

- After the operation, remove the MAIN HANDLE from the surgical lamp.
- Remove visible foreign substances using a soft brush and cloth in warm water flowing through the MAIN HANDLE for at least 1 minute.
- Place the MAIN HANDLE in the cleaning equipment and carry out the cleaning using an enzyme detergent according to the following cleaning conditions.

No.	Condition	Temperature(°C)	Minimum time(Min)
1	Pre-cleaning	20-40	1
2	Cleaning	20-40	5
3	Cleaning	20-40	2
4	Rinsing	20-40	10
5	Rinsing	20-40	10
6	Dry	50-70	90

- Remove the MAIN HANDLE from the cleaning equipment and ensure that the cleaning is complete.
- Check for foreign substances left inside and outside the MAIN HANDLE, if it is necessary, repeat the cleaning process
- Protect the cleaned and dried MAIN HANDLE from re-contamination.

## 9.3 MAIN HANDLE MANUAL CLEANING

- After the operation, remove the MAIN HANDLE from the surgical lamp.
- Immerse the MAIN HANDLE in the washing tank for at least 5 minutes in the diluted enzyme detergent. (Follow the manufacturer's instructions for use of the cleaner.)
- To prevent blood or foreign substances from drying out, soak the separated MAIN HANDLE in distilled water or tap water for 15 minutes to remove any foreign substances, then clean it with a soft brush and a lint-free cloth.
- Check inside and outside of the MAIN HANDLE for debris, and repeat cleaning if any debris remains.
- Rinse thoroughly with clean water and wipe clean with a lint-free cloth to dry.
- Protect the cleaned and dried main handle from re-contamination

**9.4 MAIN HANDLE Sterilization**

**Warning**

Do not allow foreign substances to enter the MAIN HANDLE during sterilization. MAIN HANDLE is guaranteed for up to 100 sterilization cycles when the same sterilization conditions are applied.

In the case of MAIN HANDLE, sterilization may cause wear, cracks, and discoloration. If these signs are detected, stop using the handle and replace it.

The MAIN HANDLE being stored must be kept clean and sterile before use.

If you are not wearing a sterile suit, be careful and do not touch the sterilized MAIN HANDLE.

- Sterilize the cleaned MAIN HANDLE according to the process below.
- Make sure the MAIN HANDLE is clean, wrap the HANDLE in sterile packaging (double packaging or equivalent) and seal it.
- When putting the MAIN HANDLE into the sterilizer, make sure that the hole is facing down to allow water to flow down.
- Place the MAIN HANDLE in the Pressurized Steam Sterilizer (AUTO CLAVE) and proceed with the sterilization cycle according to the sterilization conditions. Using chemical sterilization is prohibited
- Follow the sterilizer manufacturer's instructions and sterilization conditions of your country.
- When storing after sterilization, store in a tray in a sterile packaging material (room temperature storage)

<b>Sterilization Condition</b>	
<ul style="list-style-type: none"> <li>• Operating Condition</li> <li>– Temperature : 132 °C</li> <li>– Pressure : 160 ± 50kPa (1.6 ±0.5 kgf/cm<sup>2</sup>)</li> <li>– Time : 10 min</li> <li>– Dry : 16 min</li> </ul>	 <b>Autoclave</b>
<b>Sterile packaging</b>	
<ul style="list-style-type: none"> <li>• Sterile Clothing</li> <li>– KIMGUARD sterile clothing</li> <li>– Model : KC500</li> </ul>	

## 10. Maintenance

### WARNING

A electrical and mechanical check-up should be done every year. Disconnect the pendant system from the mains before any maintenance work to prevent electric shock.

### WARNING

Check points on maintenance work

- Defects of paint, Plastic fissures, Loosened parts.
- Free rotatability / limit stops and deformation of the suspension.
- The connection between LIGHTHEAD and pendant system.
- Defect of the HEAD LOCKING KEY. / The faultless function of the LIGHTHEAD.

### WARNING

In case of failure or damage, please contact your supplier.

## 11. Troubleshooting


### 11.1 LIGHTHEAD

No.	Problem	Cause	Corrective action
1	If the LED light on the LIGHTHEAD does not turn on.	Power cut	Check if supply mains are operating.
		Does not working of a backup power supply(UPS).	Please contact your supplier.
		Others	If the condition does not improve, do not attempt to repair. Please contact your supplier.
2	If there is an error in the light pattern formation.	Inappropriate distance	Check if recommended distance(1m) between the LIGHTHEAD to the bed. If the condition does not improve, do not attempt to repair. Please contact your supplier.
3	If the SAFETY COVER of the LIGHTHEAD is polluted.	Pollution	Cleanse it with designated chemicals (Alcohol, Ethanol). If the essential performances of the device (Intensity of Illumination, Color Temperature) are seriously impaired, please contact your supplier.


### 11.2 Android application(Optional)

#### 11.2.1 If the Wi-Fi connection is lost

- Check the Wi-Fi connection
- If the Wi-Fi connection is lost, please reconnect. (See section 7.3.4)
- Check the Wi-Fi Connection(See section 7.3.4)
- Reboot LIGHTHEAD and TABLET PC
- Remove DENTISLIGHT Wi-Fi and Reconnect DENTISLIGHT Wi-Fi


<TABLET Wi-Fi>


<Remove Wi-Fi>

### 11.2.2 If there is no Wi-Fi signal


- Check that the Room No. in TABLET, Wi-Fi and LIGHTHEAD controller are the same.

<TABLET Wi-Fi>

<Setting page in TABLET app>

<Setting page in LIGHTHEAD controller>


- If the Room No. is different, change the same number based on LIGHTHEAD controller and save it.
- Reboot LIGHTHEAD and TABLET
- Wi-Fi reconnection
- Check the Wi-Fi state(See section 7.3.3)

**WARNING**

If none of the troubleshooting steps above solve your problem, it may need to replace the part. Please contact your supplier or DENTIS.

## 12. Disposal




- For environment and safety of human, wastes must be recycled or separated.
  - The materials should be carefully separated.
  - The electrical boards should be submitted to an appropriate recycling proceeding
  - The cardboard box may be recycled with other paper products.
- Please contact the closest DENTIS branch or your supplier, if you have any questions about recycling of the device.

### 13. Model designation

Model	-	ARM Axis	Mount Type	Camera	Wall Controller	Battery
S300	-	2=2 Axis 3=3 Axis(SA300)	0=Head Only 1=Ceiling Type 2=Dual Ceiling Type 3=Wall Type 4=Dual Connection Type 5=Mobile Type 6=Dual ceiling(+monitor arm) 7=Triple(+monitor arm)	0=None 1= Camera	0=None	0 = None 1 = 1 EA 2 = 2 EA

**14. List of component**


 <p>CEILING VERTICAL ARM</p>	 <p>DUAL CEILING VERTICAL ARM &amp; DUAL BRACKET</p>
 <p>WALL SECOND ARM</p>	 <p>TENSION SPRING ARM</p>

 <p>MOBILE VERTICAL ARM</p>	 <p>SMPS</p> <p>+</p>  <p>SMPS MODULE</p>
 <p>S300 LIGHTHEAD</p>	 <p>WALL BRACKET</p>
 <p>CHAIR SECOND ARM</p>	 <p>CEILING SECOND ARM</p>
 <p>ARM BOLT &amp; DUAL CONNECTOR</p>	 <p>FOOT SWITCH(OPTIONAL)</p> <p>- Not applicable for export to the United States and CANADA</p>

Professional LED Surgical Luminaire for All of Surgery Application

CEILING COVER & BRACKET


MOBILE BASE & CASTER

SMPS CASE COVER



TABLET PC(OPTIONAL)

REMOTE CONTROLLER

### 15. Manufacturer and Factory


- Manufacturer, Factory : DENTIS MEDICAL DIVISION  
6, Yuram-ro, Dong-gu, Daegu, Republic of Korea

Professional LED Surgical Luminaire for All of Surgery Application

## 16. Europe Representative

- Location of EC REP  
KTR Europe GmbH  
Address : Mergenthalerallee 77, 65760 Eschborn, Germany

EL2SNP-PM080  
(V2.3) Nov/01/2024  
USR-S3-02

[www.luvis.co.kr](http://www.luvis.co.kr)

---

**DENTIS Co., Ltd. (HEADQUARTERS)** 87, Yuram-ro Dong-gu, Daegu, Republic of KOREA Tel +82-53-583-2804 Fax +82-53-583-2806  
**DENTIS MEDICAL DIVISION(FACTORY)** 6, Yuram-ro, Dong-gu, Daegu, Korea Tel. +82-70-7730-4146 Email. [luvis@luvis.co.kr](mailto:luvis@luvis.co.kr)